



Body Composition Monitor and Scale Comparison Chart



Model Number	Body Composition Monitor and Scale HBF-514C	Body Composition Monitor and Scale HBF-516B	Body Composition Monitor and Scale BCM-500
Body Weight Range	0 to 330 lbs. with an increment of .2 lbs.	0 to 330 lbs. with an increment of .2 lbs.	4.4 to 330 lbs. with an increment of .2 lbs.
Body Fat Percentage	5.0 to 60.0% with an increment of .1%	5.0 to 60.0% with an increment of .1%	5.0 to 60.0% with an increment of .1%
Body Mass Index	7.0 to 90.0 with an increment of 0.1	7.0 to 90.0 with an increment of 0.1	7.0 to 90.0 with an increment of 0.1
Skeletal Muscle Percentage	5.0 to 50.0% with an increment of .1%	5.0 to 50.0% with an increment of .1%	5.0 to 50.0% with an increment of .1%
Visceral Fat Level	30 levels with an increment of 1 level	30 levels with an increment of 1 level	30 levels with an increment of 1 level
Resting Metabolism	385 to 3999 kcal with an increment of 1 kcal	385 to 3999 kcal with an increment of 1 kcal	385 to 3999 kcal with an increment of 1 kcal
Body Age	18 to 80 year old with 1 year increments	18 to 80 year old with 1 year increments	
Memory	90-Day	180-Day	120 Total Readings in Unit Unlimited in HeartAdvisor App
Bluetooth			✓
User Profiles	4	4	4
Sensor Technology	Full Body	Full Body	Foot - to - Foot
Battery	4 AA	4 AA	4 AAA