

# OMRON Nebulizer Comparison Chart

Find out which one is right for your child. Ask your healthcare provider for a prescription to complete your purchase.



**Portable MicroAir®  
Battery-Operated  
Nebulizer**



**CompAir®  
Compressor  
Tabletop  
Nebulizer System**



**CompAir®  
Compressor  
Nebulizer with  
Kid's Accessory**

<b>NDC Number</b>	73796-0451-00	73796-0458-01	73796-8001-47
<b>Handheld</b>	✓		
<b>Tabletop</b>		✓	✓
<b>Kid-friendly design</b>			✓
<b>For children and adults</b>	✓	✓	
<b>Sound</b>	<23 dBA; virtually silent operation	≈52dBA of sound while in use, similar to the sound of an electric toothbrush*	≈45dBA of sound while in use, a little more than a quiet library*
<b>Power</b>	AA Battery or AC Adapter	AC Adapter	AC Adapter
<b>Weight</b>	4.2 oz	9.5 oz; 3.5x lighter than most tabletop compressor models	6.3 oz; our lightest tabletop compressor model
<b>360° operation for use in any position</b>	✓		
<b>Mass Mean Aerodynamic Diameter (MMAD)</b>	5 μm	3 μm	3 μm
<b>Warranty</b>	2 years	3 years	3 years
<b>Reimbursable</b>	✓	✓	✓

\*dBA comparisons. Center for Hearing and Communication. Retrieved August 19, 2020, from <https://chcheating.org/noise/common-environmental-noise-levels/>

© 2020 Omron Healthcare, Inc. All rights reserved.